

Checking Blood Glucose Levels With An Accu-Chek® Active Glucose Monitor



Insert coding chip into the side of the monitor.



Insert test strip with the orange colored square facing up. Make sure the code number on the bottle matches the code number on the machine. Pull off the tip of the lancet device.



Insert the lancet. Then, replace the tip. Select the depth number on the lancet device by turning the tip. The higher the number, the deeper the lancet goes into the skin.



Clean the area with alcohol that will be pricked. Then, push in the plunger on the lancet device. Place lancet device on finger and press button.



After the number code goes away, it will show a drop of blood on the screen. Squeeze the finger until there is a big drop of blood. Place the drop of blood on top of the orange square on the strip. After a few seconds the result will appear. Discard the strip. The machine will turn off automatically.

