

## Checking Blood Glucose Levels With A BD Logic™ Glucose Monitor

1. Insert test strip with the shiny part facing down into the top of the glucose monitor. The side with BD on it should face up. Make sure the code on the screen matches the code on the bottle. If it does not match, change the code by pressing either arrow until the numbers match. Once it finishes coding the drop of blood will show on the screen.



2. Unscrew the top of the device. Insert the lancet. Screw top back on.



3. Select the depth number on the lancet device by turning the rough area of the device. The higher the number, the deeper the lancet goes into the skin. Clean the finger to be pricked with alcohol. Pull back the end of the lancet device.



4. Place on finger and press the blue button. Squeeze the finger until there is a drop of blood. Place drop of blood on the top of the strip. After 5 seconds, it will show your

